NON-CONTACT INFRARED FOREHEAD THERMOMETER

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DEAR DR. FREI NON-CONTACT INFRARED FOREHEAD THERMOMETER OWNER

Thank you for choosing TM **Dr. Frei**® Non-Contact Infrared Forehead Thermometer Model MI-100. We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Swiss Trademark **Dr. Frei**®.

This innovative medical device uses advanced infrared (IR) technology to measure temperature instantly and accurately on the forehead or object.

Before starting to use this device please study the user's manual carefully. The user's manual offers all information you need to correctly measure your blood pressure and pulse. For all questions concerning the device please contact your local distributor or **Dr. Frei**® service center in your country.


TABLE OF CONTENTS

IMPORTANT INFORMATION ON BODY TEMPERATURE AND ITS MEASUREMENT	6
About Normal Body Temperature & Fever	6
ADVANTAGES OF NON-CONTACT INFRARED FOREHEAD THERMOMETER MODEL MI-100	7
Important information	7
BEFORE YOU BEGIN	8
Cautions and Warnings	8
Meter Overview and Display Screen	8
Select the Display Mode.....	10
Select the Temperature Unit.....	11
Replacing the Battery	11
PERFORM THE TEST	13
Using the Device As a Body Measurement Thermometer	13
Using the Device As an Object Measurement Thermometer	15
METER MEMORY	17
Recalling the Memory.....	17
Measurement Result	18
Switch OFF the Device Manually	18
SYSTEM TROUBLESHOOTING	19
MAINTENANCE	19
SYMBOL INFORMATION	20
SPECIFICATIONS	21
WARRANTY	22
APPENDIX QUICK GUIDE	102

IMPORTANT INFORMATION ON BODY TEMPERATURE AND ITS MEASUREMENT

About Normal Body Temperature & Fever

The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the Dr.Frei Non-Contact Infrared forehead thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.

Body temperature can vary from one individual/person to next. It also varies by location on the body and time of day. Below shows the statistical normal ranges from different sites. Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared. Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, overdressing or immunization. Some people may not experience fever even when they are ill. These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. Please consult your physician when you feel ill even if you do not have fever.

Table *1 Body Site Normal Temperature Range

Body Site	Normal Temperature Range
Oral	0.6°C (1°F) or more above or below 37°C (98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (armpit)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

* Body Temperature at WebMD;
website: <http://firstaid.webmd.com/body-temperature> retrieved at 2010 Jan 7.

ADVANTAGES OF NON-CONTACT INFRARED FOREHEAD THERMOMETER MODEL MI -100

Important information

This thermometer is clinically proven to produce accurate temperature measurements.

The thermometer model MI-100 is intended for the intermittent measurement and monitoring of human body temperature from forehead measurement at home.

The thermometer measures the infrared heat generated by the surface of the skin over the vessel and its surrounding tissue.

The thermometer then converts it into a temperature value.

- Automatically power off if left idle for 25 seconds and the backlight switch off in first 5 seconds.
- Memory function allows you to store 20 most recent body mode temperature results.
- Easy to read LCD with backlight in a dark environment.
- Color visible indication of alarm (red) and normal reading (green).
- Distance between device and forehead measurement within 5 cm.

NOTE:

Please be advised that the accuracy cannot be ensured when the thermometer is not clean. Check that the detecting sensor is clean before taking a measurement

A control measurement using a conventional thermometer is recommended in the following cases:

1. If the reading is surprisingly low.
2. For new-born infants up to 100 days old.
3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever.

The thermometer does not emit any infrared signal.


BEFORE YOU BEGIN


Cautions and Warnings

1. As with any thermometer, proper technique is crucial to obtaining accurate temperature readings. Please read this manual thoroughly and carefully before using.
2. Always operate the thermometer in an operating temperature range 16°C to 40°C (60.8°F to 104°F), and relative humidity less than 95%.
3. Always store the thermometer in a cool and dry place -20°C to 50°C (-4°F to 122°F) and relative humidity less than 95%.
4. Avoid direct sunlight.
5. Avoid dropping the thermometer, if it happens and you think the thermometer may be damaged, please contact customer services immediately.
6. Do not disassemble the thermometer.
7. Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.
8. This thermometer is not intended to substitute for a consultation with your physician.
9. This thermometer and the subject must remain in a stable environment for at least 30 minutes before measuring the temperature.
10. When the measured temperature falls within the temperature range of $\geq 37.5^{\circ}\text{C}$ (99.5°F) and $< 43^{\circ}\text{C}$ (109.4°F), as indicated by the red backlight on display, please consult with your physician immediately.

Meter Overview and Display Screen


Package Contents:

1. On/Scan Button
2. Display Screen
3. Guiding Light

4. Detecting Sensor
5. Battery Compartment


1. Object temperature indicator
2. Body temperature indicator
3. Low battery indicator
4. Memory mode
5. Memory set
6. Temperature unit
7. Temperature reading

Select the Display Mode

Two modes can be selected.

1. Body Mode.


Press On/Scan button for less than 2 seconds to turn on the thermometer first and entry into body mode. When ° C or ° F blinks, the device is standby for measurement.

This mode is used to measure the forehead temperature.

2. Object Mode.

Press On/Scan button for 2~5 seconds to turn on the thermometer first and entry into object mode. When ° C or ° F blinks, the device is standby for measurement.


This mode is used to measure the object temperature.

Select the Temperature Unit

This meter provides two measurement units used for indicating the body/object temperature, °C or °F, for your preferred selection.

When device is off, press On/Scan button for 8 seconds (the memory recall interface will display before entering into unit selection operation) to set up the unit, as °C or °F is blinking, press the On/Scan button to switch between unit C and F. The device is standby for measurement


The device will switch off since no button pressed within 5 seconds.

Replacing the Battery

The thermometer comes with one 3V CR2032 lithium battery. The meter will display «» to alert you when the meter power is getting low, please follow the steps below to replace new battery immediately.


1. Press the buckle on battery cover and lift up to Remove the cover.


2. Remove the old battery with battery hook by lifting the battery slightly upward and replace with one new 3V CR2032 lithium battery.

3. Close the battery cover.

NOTE:

1. Although the thermometer works when «» appears, we still recommend that you change the battery to obtain an accurate result.
2. Remove the battery if stored for a long period of time.
3. The battery should be kept away out of child's reach.
4. The battery shall be properly disposed according to your local environmental regulations.
5. Follow the national requirement to dispose the unit.

PERFORM THE TEST

Using the Device As a Body Measurement Thermometer

1. Press On/Scan button for less than 2 seconds to turn on the thermometer first and entry into body mode.


The full display will be showed on the meter for 2 seconds.

2. The last memory displays for 2 seconds.

3. When ° C or ° F blinks, the device is standby for measurement. The temperature unit flashes.

4. Move the detecting sensor close to the forehead. Make sure the detecting sensor is flat and close to the forehead, not at an angle. Perform a forehead measurement with a distance within 5 cm.

Press On/Scan button for up to 3 seconds, the meter will take measurement.

5. Read the result.

The measurement result will be done once releasing the button. The reading is shown together with green backlight and three short beeps.


NOTE:

- As the forehead measurement temperature is likely to be affected by sweat, oil and the surroundings, the reading shall be taken as a reference only.
- Babies' skin reacts very quickly to the ambient temperature. Therefore, do not take their temperature with the Dr.Frei Non-Contact Infrared forehead thermometer during or after breastfeeding because the skin temperature may then be lower than the internal body temperature.
- If the reading is $< 32^{\circ}\text{C}$ (89.6°F), the display will show "L" together with red backlight and three short beeps.
- If the reading is $\geq 32^{\circ}\text{C}$ (89.6°F) and $< 37.5^{\circ}\text{C}$ (99.5°F), the display will show together with green backlight and three short beeps.
- If the reading is $\geq 37.5^{\circ}\text{C}$ (99.5°F) and $< 43^{\circ}\text{C}$ (109.4°F), the display will show together with red backlight and ten short beeps.
- If the reading is $\geq 43^{\circ}\text{C}$ (109.4°F), the display will show "H" together with red backlight and ten short beeps.
- The interval for the next measurement should be 3 seconds.
- Device is standby for another measurement with last readings stays on the display, $^{\circ}\text{C}$ or $^{\circ}\text{F}$ is blinking. The device cannot do another measurement during previous measurement shows on the display with backlight that lasts for 5 sec.
- The thermometer will automatically turn off if left idle for 25 seconds and the backlight switch off in first 5 seconds.

Using the Device As an Object Measurement Thermometer

1. Press On/Scan button for 2~5 seconds to turn on the thermometer first and entry into object mode.

The object mode icon will appear with a short beep during this operation. The full display will be showed on the meter for 2 seconds.

2. The last memory displays for 2 seconds.

3. When ° C or ° F blinks, the device is standby for measurement. The temperature unit flashes.

4. Move the detecting sensor close to the object.

Make sure the detecting sensor is flat and close to the object, not at an angle. Perform an object measurement with a distance within 5 cm.

Press On/Scan button for up to 5 seconds, the meter will take measurement.

5. Read the result.

The measurement result will be done once releasing the button. The reading is shown together with green backlight and three short beeps.

NOTE:

- If the reading is $< 0^{\circ}\text{C}$ (32°F), the display will show "L" together with red backlight and three short beeps.
- If the reading is $\geq 0^{\circ}\text{C}$ (32°F) and $< 100^{\circ}\text{C}$ (212°F), the display will show together with green backlight and three short beeps.
- If the reading is $\geq 100^{\circ}\text{C}$ (212°F), the display will show "H" together with red backlight and ten short beeps.
- The interval for the next measurement should be 3 seconds.
- Device is standby for another measurement with last readings stays on the display, $^{\circ}\text{C}$ or $^{\circ}\text{F}$ is blinking. The device cannot do another measurement during previous measurement shows on the display with backlight that lasts for 5 sec.
- The thermometer will automatically turn off if left idle for 25 seconds and the backlight switch off in first 5 seconds.

METER MEMORY

Recalling the Memory

Your thermometer stores 20 most recent body temperature readings.

1. Be sure the thermometer is OFF before recalling this memory.
2. Press On/Scan button for 5~8 seconds to enter the memory mode.

Each time you press the On/Scan button, a result will be displayed in the order (latest result shown first), together with «» and number (from 1 to 20).

Backlight in green or red will appear according to the memory reading.

When the memory is full, the oldest result is deleted as the new one added. When the last record displayed in the Display, press Memory Button again to return the first record.

If there is no readings existed, «» will be shown on the display while «» is blinking.

3. Exit the memory.

The device will switch off since no button pressed within 5 seconds.

NOTE:

Only the body temperature records can be read in this memory mode.

Measurement Result



Message	What it means
	Measured temperature falls outside the displayed temperature range of 32°C to 43°C (89.6°F to 109.4°F), is shown on the display with red backlight:
	<ul style="list-style-type: none">• Hi-temperature $\geq 43^{\circ}\text{C}$ (109.4°F)• Lo-temperature $< 32^{\circ}\text{C}$ (89.6°F)
	Measured temperature falls in normal reading, within the range of 32°C to 37.5°C (89.6°F to 99.5°F), is shown on the display with green backlight.
	Measured temperature falls within the range of $\geq 37.5^{\circ}\text{C}$ (99.5°F) and $< 43^{\circ}\text{C}$ (109.4°F), is shown on the display with red backlight.

Switch off the Device manually

After the measurement is completed, Press On/Scan button to turn off the thermometer while the backlight is on, either using as a body measurement thermometer or an object measurement thermometer.

SYSTEM TROUBLESHOOTING

Message	What it means	What to do
	Appear when environmental temperature is below system operation range, or the object temperature is below displayed temperature range.	Put the thermometer under operating temperature range of 16°C to 40°C (60.8°F to 104°F), or under object temperature range of 0°C to 100°C (32°F to 212°F)
	Appear when environmental temperature is above system operation range, or the object temperature is above displayed temperature range.	Put the thermometer under operating temperature range of 16°C to 40°C (60.8°F to 104°F), or under object temperature range of 0°C to 100°C (32°F to 212°F)
	Battery is low and «  » appears on display.	Please replace battery as soon as possible.
	If the thermometer is broken or affected by EMF.	Please contact the customer service line.

MAINTENANCE

Cleaning and Care

- The detecting sensor is not waterproof. Please wipe with a clean and dry cotton swab to clean the detecting sensor on the inside.

- The body of the thermometer is not water-resistant. Never put the thermometer under a running tap or submerge it into water. Use a soft and dry cloth to clean it. Do not use abrasive cleaners.
- Store the thermometer in a cool and dry location. Free from dust and away from direct sunlight.

SYMBOL INFORMATION

Symbol	Referent
	Consult instructions for use
	Manufacturer
SN	Serial number
	Caution, consult accompanying documents
CE ₀₁₉₇	Mark CE
	Authorized representative in the European Community
	Type BF Equipment
	Collection for electrical and Electronic equipment

SPECIFICATIONS

Model No	MI-100
Dimensions & Weight	100 x 30.5 x 18 mm, 38 g (include batteries)
Power Source	One 3V CR2032 lithium battery
Displayed temperature range	Forehead: 32°C to 43°C (89.6°F to 109.4°F) Object: 0°C to 100°C (32°F to 212°F)
Display resolution	0.1 °F/0.1°C
Accuracy	For body mode: ±0.2°C (±0.36°F) from 36.0°C (96.8°F) to 39.0°C (102.2°F) ±0.3°C (±0.54°F) from 32.0°C (89.6°F) to 35.9°C (96.6°F) from 39.1°C (102.4°F) to 43.0°C (109.4°F) For object mode: ± 1°C (±2°F) from 0°C (32°F) to 100°C (212°F)
Temperature unit	°C or °F
Operating temperature range	16°C to 40°C (60.8°F to 104°F)
Operating humidity	95% RH or less
Storage temperature range	-20°C to 50°C (-4F to 122°F)
Storage humidity	95% RH or less
Memory capacity	20 measurements
Standards	Device Corresponds to the requirements of the standard for infrared thermometers ASTM E1965-98 EN 12470-5 EN 60601-1-2:2007/AC:2010 EN 60601-1-4:1996 EN 60601-1-6:2010

	<p>Electromagnetic Compatibility: Device fulfills the stipulations of the standard EN 60601-1-2</p> <p>The stipulations of EU-Directive 93/42/EEC for Medical Devices Class IIa have been fulfilled.</p>
--	--

* The specifications and the outer design are subject to change for improvement reasons without previous notice.

Batteries and electronic instruments must be disposed Of in accordance with the locally applicable regulations, Not with domestic waste.

WARRANTY

This product is under warranty for 2 years after date of purchase. The warranty is only valid upon presentation of the warranty card which was correctly filled in and sealed by the seller (or by the service center representative) with the date of the purchase (repair), the check, instruction manual and safe packaging of the device.

- The warranty does not apply to wear out of battery and packaging box.
- The warranty does not apply to damage caused by improper handling, accidents, not following the operating instructions or self-maintained alterations made to the device.
- The device cannot be returned if it is without original packaging, with mechanical damages (scratches, stains and so on) or without any component.
- The device which was used and needs repair cannot be replaced with a new one.

In case there was violation of operating rules during warranty period, repair will be done at the expense of the consumer.

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Manufactured by the order of: MedHouse Swiss GmbH
■ Biocare Asia Corporation, No.260 Mayun Road, New District Suzhou
215129 Jiangsu, China.   Shanghai International Trading Corp. GmbH